



smith&nephew

EVOS[◇] SMALL

Plating System



Supporting healthcare professionals

2.7mm		4.0mm		3.5mm		4.7mm	
Cortex	Locking	Fully Threaded Osteopenia	Partially Threaded Osteopenia	Cortex	Locking	Fully Threaded Osteopenia	Partially Threaded Osteopenia



Thread diameter	2.7mm	2.7mm	4.0mm	4.0mm	3.5mm	3.5mm	4.7mm	4.7mm
Head diameter	4.5mm	4.3mm	4.5mm	4.5mm	5.6mm	5.4mm	5.6mm	5.6mm
Core diameter	2.0mm	2.0mm	2.0mm	2.0mm	2.5mm	2.5mm	2.5mm	2.5mm
Thread pitch	1.0mm	1.0mm	1.75mm	1.75mm	1.25mm	1.25mm	1.75mm	1.75mm
Driver	T8	T8	T8	T8	2.5mm Hex	2.5mm Hex	2.5mm Hex	2.5mm Hex
Screw lengths	6–22mm (1mm Increments)	6–22mm (1mm Increments)	10–50mm (2mm Increments)	26–50mm (2mm Increments)	6–20mm (1mm Increments)	8–20mm (1mm Increments)	10–50mm (2mm Increments)	26–50mm (2mm Increments)
	24–50mm (2mm Increments)	24–50mm (2mm Increments)	55–80mm (5mm increments)	55–80mm (5mm increments)	22–50mm (2mm Increments)	22–50mm (2mm Increments)	55–110mm* (5mm Increments)	55–110mm* (5mm Increments)
	55–80mm (5mm increments)	55–80mm (5mm increments)			55–110mm* (5mm Increments)	55–110mm* (5mm Increments)		
Thread length	-	-	-	40% of screw length	-	-	-	40% of screw length

Mediala distala Tibia

Partial Articular 2.7/3.5mm Medial Distal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72453903	72454003	3	64	
72453906	72454006	6	97	
72453909	72454009	9	130	



Standard 2.7/3.5mm Medial Distal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72463909	72464009	9	130	1
72463912	72464012	12	162	1
72463915	72464015	15	195	1
72463918	72464018	18	228	
72463921	72464021	21	261	



Instrument: EVOS small

Kortikal & vinkelstabil skruv 2.7 & 3.5mm

Spongiös skruv: 4,0 & 4.7 mm

Anterior Distal Tibia

Partial Articular 2.7/3.5mm Anterior Distal Tibia				
Art nr	Hål	Längd mm		Har hemma
72454103	3	74		1
72454106	6	107		1
72454109	9	140		



Instrument: EVOS small

Kortikal & vinkelstabil skruv 2.7 & 3.5mm

Spongiös skruv: 4,0 & 4.7 mm

Anterolateral Distal Tibia

Partial Articular 2.7/3.5mm Anterolateral Distal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72454306	72454406	6	87	1
72454309	72454409	9	120	



Standard 2.7/3.5mm Anterolateral Distal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72464308	72464408	8	120	1
72464311	72464411	11	153	1
72464314	72464414	14	186	1
72464317	72464417	17	219	
72464320	72464420	20	252	



Instrument: EVOS small

Kortikal & vinkelstabil skruv 2.7 & 3.5mm

Spongiös skruv: 4,0 & 4.7 mm

Posterior Distal Tibia

Partial Articular 3.5mm Posterior Distal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72454503	72454603	3	63	1
72454506	72454606	6	96	1



Standard 3.5mm Posterior Distal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72464506	72464606	6	98	
72464509	72464609	9	131	
72464512	72464612	12	163	
72464515	72464615	15	196	



Instrument: EVOS small

Kortikal & vinkelstabil skruv 3.5mm

Spongiös skruv: 4.7 mm

Lateral Proximal Tibia

Partial Articular 3.5mm Lateral Proximal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72453104	72453204	4	70	
72453106	72453206	6	91	



Standard 3.5mm Lateral Proximal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72463104	72463204	4	70	1
72463106	72463206	6	91	1
72463108	72463208	8	113	1
72463110	72463210	10	134	1
72463113	72463213	13	167	1
72463116	72463216	16	200	
72463118	72463218	18	221	
72463120	72463220	20	243	
72463122	72463222	22	265	
72463124	72463224	24	287	




Instrument: EVOS small

Kortikal & vinkelstabil skruv 3.5mm

Spongiös skruv: 4.7 mm

Medial Proximal Tibia

Partial Articular 3.5mm Medial Proximal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72453704	72453804	4	75	
72453708	72453808	8	117	



Standard 3.5mm Medial Proximal Tibia				
Vänster	Höger	Hål	Längd mm	Har hemma
72463708	72463808	8	117	1
72463710	72463810	10	138	1
72463713	72463813	13	170	
72463716	72463816	16	201	



Instrument: EVOS small

Kortikal & vinkelstabil skruv 3.5mm

Spongiös skruv: 4.7 mm

Posteromedial proximal Tibia

Partial Articular 3.5mm Posteromedial Proximal Tibia "T"				
Vänster	Höger	Hål	Längd mm	Har hemma
72453304	72453404	4	71	1
72453307	72453407	7	103	1



Standard 3.5mm Posteromedial Proximal Tibia "T"				
Vänster	Höger	Hål	Längd mm	Har hemma
72463307	72463407	7	104	1
72463310	72463410	10	137	1
72463313	72463413	13	170	1



Partial Articular 3.5mm Posteromedial Proximal Tibia "I"				
Vänster	Höger	Hål	Längd mm	Har hemma
72453505	72453605	5	78	
72453508	72453608	8	111	



Instrument: EVOS small

Kortikal & vinkelstabil skruv 3.5mm

Spongiös skruv: 4.7 mm

Lateral Distal Fibia

2.7mm Lateral Distal Fibula				
Vänster	Höger	Hål	Längd mm	Har hemma
72464705	72464805	5	61	
72464708	72464808	8	82	
72464711	72464811	11	103	



2.7/3.5mm Lateral Distal Fibula				
Vänster	Höger	Hål	Längd mm	Har hemma
72465103	72465203	3	59	
72465105	72465205	5	81	2
72465107	72465207	7	103	3
72465109	72465209	9	125	3
72465111	72465211	11	147	2
72465113	72465213	13	169	1
72465116	72465216	16	202	



3.5mm Lateral Distal Fibula				
Vänster	Höger	Hål	Längd mm	Har hemma
72464903	72465003	3	59	
72464905	72465005	5	81	2
72464907	72465007	7	103	2
72464909	72465009	9	125	1
72464911	72465011	11	147	
72464913	72465013	13	169	1
72464916	72465016	16	202	



Instrument: EVOS small

Kortikal & vinkelstabil skruv 2.7 och 3.5mm

Spongiös skruv: 4,0 & 4.7 mm

Posterolateral Distal Fibula

3.5mm Posterolateral Distal Fibula Anti Glide				
Vänster	Höger	Hål	Längd mm	Har hemma
72455305	72455405	5	58	1
72455306	72455406	6	69	1
72455307	72455407	7	80	



2.7/3.5mm Posterolateral Distal Fibula				
Vänster	Höger	Hål	Längd mm	Har hemma
72465305	72465405	5	93	1
72465307	72465407	7	115	1
72465309	72465409	9	137	1
72465311	72465411	11	159	
72465314	72465414	14	192	



Instrument: EVOS small

Kortikal & vinkelstabil skruv 2.7 och 3.5mm

Spongiös skruv: 4,0 & 4.7 mm

Raka plattor 2.7mm

2.7mm Kompression platta			
Art nr	Hål	Längd mm	Har hemma
72440404	4	33	
72440406	6	50	1
72440408	8	67	1
72440410	10	84	1
72440415	15	127	
72440418	18	153	



2.7mm Vinkelstabil Kompressions platta			
Art nr	Hål	Längd mm	Har hemma
72440204	4	31	
72440206	6	46	1
72440208	8	61	1
72440210	10	76	1
72440215	15	113	
72440218	18	136	



2.7mm Recon platta			
Art nr	Hål	Längd mm	Har hemma
72440304	4	33	
72440306	6	49	
72440308	8	65	
72440310	10	81	
72440315	15	121	
72440318	18	145	



2.7mm Vinkelstabil Recon Platta			
Art nr	Hål	Längd mm	Har hemma
72440104	4	32	
72440106	6	48	
72440108	8	64	
72440110	10	80	
72440115	15	120	
72440118	18	144	



Raka Plattor 3.5mm

3.5mm Kompression Platta			
Art nr	Hål	Längd mm	Har hemma
72441004	4	52	
72441006	6	77	5
72441007	7	90	6
72441008	8	102	5
72441010	10	127	2
72441012	12	152	2
72441014	14	177	
72441016	16	202	
72441018	18	227	
72441020	20	252	



3.5mm Vinkelstabil Kompression Platta			
Art nr	Hål	Längd mm	Har hemma
72440704	4	47	1
72440706	6	70	1
72440707	7	81	1
72440708	8	93	1
72440710	10	116	1
72440712	12	139	1
72440714	14	162	1
72440716	16	185	
72440718	18	208	
72440720	20	231	



3.5mm Recon Plate			
Art nr	Hål	Längd mm	Har hemma
72440904	4	44	
72440906	6	66	1
72440908	8	88	1
72440910	10	110	1
72440912	12	132	
72440914	14	154	
72440916	16	176	
72440918	18	198	
72440920	20	220	
72440922	22	242	



Raka Plattor 3.5mm forts.

3.5mm Vinkelstabil Recon platta			
Art nr	Hål	Längd	Har hemma
72440604	4	44	
72440606	6	66	
72440608	8	88	
72440610	10	110	
72440612	12	132	
72440614	14	154	
72440616	16	176	
72440618	18	198	
72440620	20	220	
72440622	22	242	



3.5mm 1/3 rd Tubular platta			
Art nr	Hål	Längd mm	Har hemma
72440802	2	22	2
72440804	4	46	2
72440806	6	70	9
72440807	7	82	6
72440808	8	94	6
72440810	10	118	3
72440812	12	142	1



3.5mm Vinkelstabil 1/3 rd Tubular platta			
Art nr	Hål	Längd mm	Har hemma
72440502	2	22	
72440504	4	46	
72440506	6	70	1
72440507	7	82	1
72440508	8	94	1
72440510	10	118	1
72440512	12	142	



2.7mm Flex Plate			
Art nr	Hål	Längd mm	Har hemma
72442720	3 Hole Head 5H	28	2
72442722	3 Hole Head 9H	71	



2.7mm Flex Plate			
Art nr	Hål	Längd mm	Har hemma
72442724	4	28	
72442726	6	42	2
72442728	8	56	2
72442760	10	70	
72442762	20	140	



2.7mm Strength Plate			
Art nr	Hål	Längd mm	Har hemma
72442764	4	28	
72442766	6	42	
72442768	8	56	
72442770	10	70	
72442772	20	140	



2.7mm Straight Tine Plate			
Art nr	Hål	Längd mm	Har hemma
72442774	5	38	
72442767	7	52	



2.7mm Y-Plate			
Art nr	Hål	Längd	Har hemma
72442778	5	55	
72442780	7	69	



2.7mm Y Tine-Plate			
Art nr	Hål	Längd	Har hemma
72442782	5	52	
72442784	7	66	



2.7mm Triangle Plate			
Art nr	Hål	Längd	Har hemma
72442786	6	54	
72442788	8	68	

