


Cannulated Screw System

- Self-Tapping conical headless screws
- Diameters from 1.7mm to 7.3mm
- 1mm length increments available




Cannulated Screw System


Wire Guide/Cannula System enables easy, quick screw insertion.




Insert Wire



Drill










Insert screw



Final fixation

Shown at actual size:



	SCREWS / LENGTHS*	WIRE GUIDES	DRILL CANNULA	WIRES	DRILL BITS	COUNTERSINKS
	1.7 C17xx 8–14mm*	WGUIDE-1.7	n/a	WIRE-0.7/080	[self-drilling]	n/a
	2.3 C23xx 10–28mm*	WGUIDE-2.3	CANNULA-2.3	WIRE-0.8/120	DRILL-1.6CL	CSINK-2.3
	3.0 C30xx 14–32mm*	WGUIDE-3.0	CANNULA-3.0	WIRE-1.1/120	DRILL-2.1CL	CSINK-3.0
	3.5 C35xx 14–45mm*	WGUIDE-3.5	CANNULA-3.5	WIRE-1.1/120	DRILL-2.4CL	CSINK-3.5
	4.5 C45xx 26–44mm 44–64mm	JD700	GM730	WIRE-1.6/230	DRILL-3.0C	CSINK-4.5
	7.3 C73xx 68–112mm** CL73xx	JD700	GM730	WIRE-2.0/230	DRILL-4.7C	CSINK-7.3

*1mm Increments available

**4mm Increments available
16mm & 32mm distal thread lengths available.

Small Headless Screws (1.7, 2.3, 3.0, 3.5mm)
Large Headless Screws (4.5, 7.3mm)